


## Timeline 2 - Allogenic Rehabilitation of Iatrogenically Degraded Gastrointestinal Flora



**(A)** This is a continuation from step 4 in Timeline 1. The individual components of the intestinal flora are cultured in various media.

**(B)** The cultured intestinal flora components are mixed together in the correct proportions. The resultant mixture may be centrifuged to concentrate it.

**(C)** The concentrated cultured intestinal flora is transferred to Lyoguard trays for freeze drying in large quantities.

**(D)** The freeze dried product can be packaged in capsules for use by patients.